

Flat Bottom Dishes



				Estimated Weight (g)	Estimated Weight (g)	Estimated Weight (g)	Estimated Weight (g)	Estimated Weight (g)
Code	Capacity (ml)	Diameter (mm)	Depth (mm)	Std. Form	Rim Only (reinforced)	Base Only (reinforced)	Rim & Base (reinforced)	Std. Lid Only
GD20	20	41	20	6	7	8	9	6
GD25	25	42	23	8	9	10	12	10
GD35	35	49	25	12	13	14	16	13
GD40	40	51	26	15	17	18	20	14
GD50	50	53	29	17	19	21	25	15
GD60	60	58	30	0	21	23	27	19
GD75	75	65	30	25	27	30	33	22
GD100	100	70	33	33	35	37	40	25
GD125	125	74	34	42	44	45	48	28
GD150	150	77	41	50	53	55	60	33

				Estimated Weight (g)	Estimated Weight (g)	Estimated Weight (g)	Estimated Weight (g)	Estimated Weight (g)
GD175	175	84	40	55	59	64	68	41
GD200	200	90	43	67	71	75	80	44
GD250	250	94	42	80	82	83	88	51
GD300	300	110	45	94	98	100	105	55
GD400	400	115	48	133	137	140	148	58
GD500	500	124	50	170	179	185	200	61
GD750	750	124	75	270	290	310	350	61