


Low Profile Crucibles



*Reinforcing of the rim gives the crucible a more robust character enabling regular handling and reducing rim deformation.

Available Capacity: 8-125ml (larger sizes made to order)

Options: Lids, reinforced rims and/or bases

					Estimated Weight (g)	Estimated Weight (g)	Estimated Weight (g)	Estimated Weight (g)
Code	Capacity (ml)	Diameter (mm)	Depth (mm)	Std. Form	Rim Only (reinforced)	Base Only (reinforced)	Rim & Base (reinforced)	Std. Lid Only
GCL8	8	25	21	6	7	7	8	2
GCL10	10	27	23	7	7.5	8	9	3
GCL15	15	30	27	11	12	13	14	4
GCL20	20	33	28	15	16	17	18	5
GCL25	25	36	30	19	20	21	23	5
GCL30	30	40	32	22	24	26	27	7
GCL35	35	42	33	25	26	27	32	7

					Estimated Weight (g)	Estimated Weight (g)	Estimated Weight (g)	Estimated Weight (g)
GCL40	40	43	34	30	32	34	36	7
GCL50	50	46	37	38	40	42	45	8
GCL60	60	48	39	46	48	50	54	9
GCL70	70	52	41	50	52	54	63	12
GCL75	75	52	44	50	53	58	67	13
GCL80	80	52	47	51	54	63	72	14
GCL90	90	57	48	52	59	71	81	15
GCL100	100	57	51	59	68	73	90	17
GCL110	10	57	54	68	71	81	95	19
GCL120	120	57	57	84	88	90	98	22
GCL125	125	58	59	97	101	105	114	25